




**Discovery and targeted phosphoproteomics
in under a week:
Phosphorylation analysis using
data-independent acquisition (DIA)**

Benjamin L. Parker










**NHMRC Early Career Fellow
Diabetes & Obesity Program
Biological Mass Spectrometry Unit Manager
David James Laboratory
Garvan Institute of Medical Research, Sydney, Australia**




Case study: Insulin stimulated phosphoproteome Humphrey *et. al. Cell Metab.* 2013



Data-Dependent Acquisition (DDA)




Data-Dependent Acquisition (DDA)



- Precursor ions trigger an isolation in a 1-3 m/z window and are subjected to MS/MS.
- *Simplified* MS/MS spectra for characterization of the isolated analyte(s).
- Many precursor ions are not subjected to MS/MS and cannot be characterized = missing values between replicates

Data-Independent Acquisition (DIA)



- Precursor ions are co-isolated without *priori* knowledge

- the entire mass range (all ion fragmentation: AIF)


Masselon, *Anal. Chem* (2000); Purvine, *Proteomics* (2003); Plumb, *Rapid Commun. Mass Spectrum*. (2006); Geiger, *Mol. Cell. Proteomics* (2010)

- 1-30 *m/z* (SWATH)

Venable, *Nat. Methods* (2004); Panchaud, *Anal. Chem.* (2009); Gillet, *Mol. Cell. Proteomics* (2012); Egertson, *Nat. Methods* (2013)


- Every precursor ion is subjected to MS/MS across the desired mass range = improved repeatability between replicates

1- Targeted quantification of AKT substrate phosphorylation using SID-DIA



15 synthetic heavy labelled phosphopeptides

1- Targeted quantification of AKT substrates using SID-DIA



260 µg protein digest Lys-C/trypsin & desalt

Peptides resuspended in TiO₂ Loading Buffer containing a normalized amount of mixed heavy labelled phosphopeptides. TiO₂ enrichment

Single-shot nanoUHPLC over 120 min DIA acquisition on a Q-Exactive: 450-850 m/z; 25 m/z isolation; loop count 8 with intermittent MS1 scans. Analysis in Skyline.




Fractionated Accucore Amide-HILIC over 40 min; 8 fractions. Each fractionation DDA acquisition. Searched with Sequest HT/Percolator in Proteome Discoverer v1.4 and spectral library created in Skyline.


1- Targeted quantification of AKT substrate phosphorylation using SID-DIA

Targeted Data Extraction *Gillet, Mol. Cell. Proteomics (2012)*

DDA Spectral Library (Heavy)




DIA (Heavy)




idotp = 0.99
dotp = 0.92

DIA (Light)



idotp = 0.98
dotp = 0.93

DIA (L/H)




Targeted Data Extraction Identification Criteria

- 1- Retention time matching of spectral library to experimentally observed within 5% of the gradient length (i.e a 5 min window).
- 2- Co-elution of endogenous (L) and heavy (H) phosphopeptides.
- 3- Dot-product matching of theoretical to experimentally observed precursor-ion isotope distribution > 0.95 .
- 4- Dot-product matching of spectral library to experimentally observed product-ion distribution > 0.90 .
- 5- Matching peak shape for precursor and product ions from both L and H phosphopeptides.

13 out of 15 phosphopeptides were successfully identified
11 out of 13 phosphopeptides had L:H within 10-fold




Quantification strategy and reproducibility

Technical coefficient of variation (CV) of the entire protocol (insulin stimulated group) using precursor-ion XICs (MS1), product-ion XICs (MS2) and summed precursor- and product-ion XICs (MS1+MS2) with and without SID normalisation.






All 13 phosphopeptides were quantified with <20% CV regardless of the quantification strategy.

1- Targeted quantification of AKT substrate phosphorylation using SID-DIA



Potent inhibition (>80%) of AKT S473/4 was observed with AKTi, PI3Ki and mTOR while PDKi showed partial attenuation on AKT1 but not AKT2. SGKi had no effect on AKT phosphorylation.



Jacinto, E., et al., Cell (2006)


2- Targeted quantification of predicted AKT and mTOR substrate phosphorylation using label-free DIA

60 top predicted AKT and 60 top predicted mTOR substrates were selected from a recent large-scale SILAC experiment (Humphrey, *et. al., Cell Metab.* (2013). Machine learning:


Spectral library containing 10,285 phosphosites and 99% low response to AKT and PI3K/mTOR phosphopeptides with 3- temporal profile.

= Pareto Rankings
Targeted data extraction:


- 45 out of 98 phosphoproteins met
 - 42 out of 45 phosphoproteins met
- DIA criteria.
with <20% CV



Targeted quantification of predicted AKT and mTOR substrate phosphorylation using label-free DIA




GAB2 (GRB2-associated-binding protein 2) – docking protein with PH domain, PxxP motifs for SH2/SH3 adaptor activity and binds to numerous signalling molecules in PI3K/AKT pathway.





Humphrey, *et. al.*, *Cell Metab.* (2013)

AKT *In vitro* kinase analysis of GAB2 by DIA




In vitro kinase analysis by DIA: AKT phosphorylates GAB2 at S210 and T391




4- 14-3-3 binding of GAB2: insulin versus EGF

Phosphorylation of S210 and T391 recruits 14-3-3 and negatively regulates EGF signalling
(Brummer, *et al.*, *EMBO J* (2008))



SID-DIA-MS



Conclusions and future work

Phosphoproteomic analysis using DIA is a viable strategy for targeted analysis with and without SID normalization.

Future work:

- further work with false-discovery rates for identification**
- further work with normalization methods**
- quantification of total protein abundances**

Acknowledgements

Prof. David James
Dr. Guang Yang
Dr. Xiuquan Ma
Dr. Rima Chaudhuri
Dr. Sean Humphrey
Dr. Daniel Fazakerley
Dr. Jacqueline Stoeckli
and The James Lab

Dr. Rosa Viner
Dr. Scott Peterman
Mr. Mark O'Brien
Dr. Steve Binos
(Thermo Fisher)
Dr. Jarrett Egertson
Mr. Brendan MacLean
(MacCoss Lab)



N | H | M | R | C

ThermoFisher
SCIENTIFIC